

SIEMENS MAGNETOM Avanto syngo MR B17

\\USER\MGH\Diffusion\DT\DIFFUSION

TA: 2:00 PAT: 2 Voxel size: 1.7x1.7x5.0 mm Rel. SNR: 1.00 USER: ep2d_diff_MGH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	On
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	23
Dist. factor	20 %
Position	R3.0 A3.0 H0.0
Orientation	T > C-12.5
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	220 mm
FoV phase	100.0 %
Slice thickness	5.00 mm
TR	3120 ms
TE	88 ms
Averages	1
Concatenations	1
Filter	Raw filter, Prescan Normalize
Coil elements	HE1-4

Contrast

MTC	Off
Magn. preparation	None
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Matrix Coil Mode	Auto (Triple)
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	On
Raw filter	On
Intensity	Weak
Slope	25
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
------------------	-------------

Series Interleaved

Special sat.	None
Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System

Body	Off
HE2	On
HE4	On
HE1	On
HE3	On

Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain Atlas
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R3.0 A3.0 H0.0
Orientation	T > C-12.5
Rotation	0.00 deg
R >> L	220 mm
A >> P	220 mm
F >> H	137 mm

Physio

1st Signal/Mode	None
Resp. control	Off

Diff

Diffusion mode	MDDW
Diff. weightings	2
b-value 1	0 s/mm ²
b-value 2	1000 s/mm ²
Mosaic	On
Noise level	40
Diff. directions	30

Sequence

Introduction	On
Bandwidth	1396 Hz/Px
Free echo spacing	Off
Echo spacing	0.8 ms
EPI factor	128
RF pulse type	Normal
Gradient mode	Fast*
Sequence Mode	Product
Diff Grad Table	Single
Direction Scheme	Single
Dummy Scans	3
T2 Weighted Images	5

SIEMENS MAGNETOM Avanto syngo MR B17

FFT Scale Factor	1.00
Diff Grad Mode	XYZ

SIEMENS MAGNETOM Avanto syngo MR B17

\\USER\MGH\Diffusion\DT\DIFFUSION_HighRes

TA: 10:02 PAT: 2 Voxel size: 2.0x2.0x2.0 mm Rel. SNR: 1.00 USER: ep2d_diff_MGH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	On
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	64
Dist. factor	0 %
Position	R3.0 A3.0 H0.0
Orientation	T > C-12.5
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	8230 ms
TE	81 ms
Averages	1
Concatenations	1
Filter	Raw filter, Prescan Normalize
Coil elements	HE1-4

Contrast

MTC	Off
Magn. preparation	None
Fat suppr.	Fat sat.

Averaging mode	Long term
Reconstruction	Magnitude
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Matrix Coil Mode	Auto (Triple)
Reference scan mode	Separate

Distortion Corr.	Off
Prescan Normalize	On
Raw filter	On
Intensity	Weak
Slope	25
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
------------------	-------------

Series

Special sat.	None

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System

Body	Off
HE2	On
HE4	On
HE1	On
HE3	On

Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain Atlas
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R3.0 A3.0 H0.0
Orientation	T > C-12.5
Rotation	0.00 deg
R >> L	256 mm
A >> P	256 mm
F >> H	128 mm

Physio

1st Signal/Mode	None

Resp. control	Off

Diff

Diffusion mode	MDDW
Diff. weightings	2
b-value 1	0 s/mm ²
b-value 2	700 s/mm ²
Mosaic	On
Noise level	40
Diff. directions	60

Sequence

Introduction	On
Bandwidth	1396 Hz/Px
Free echo spacing	Off
Echo spacing	0.8 ms

EPI factor	128
RF pulse type	Normal
Gradient mode	Fast*

Sequence Mode	Product
Diff Grad Table	Single
Direction Scheme	Single
Dummy Scans	3
T2 Weighted Images	10

SIEMENS MAGNETOM Avanto syngo MR B17

FFT Scale Factor	1.00
Diff Grad Mode	XYZ

Table of contents

\\USER	MGH	Diffusion	DTI	DIFFUSION
				DIFFUSION_HighRes