

Note: To ensure correct and safe use of the products, make sure that you read the Operating Instructions on the DVD carefully.



BRAINVISION professional **BRAINAMP Set Up** ✓

1

Installing the software: Insert the Brain Products Application Suite DVD and install the Recorder software as described in the User Manual. The User Manual can be found on the Application Suite DVD.


Important: Leave the DVD in the drive while installing the hardware.



2

Installing the dongle: Plug in the USB dongle. The dongle driver is installed automatically. When installation is complete, the dongle will light up.

Important: Your computer must be equipped with at least two USB2 ports for the software license dongle and the USB2 Adapter.




3

Installing the USB2 Adapter: Use the USB cable to connect the USB2 Adapter to your computer. Windows will detect the new hardware and install it automatically.

4

Installing the amplifier: Connect the amplifier to the USB2 Adapter using the fiber optic cable (FOC5 or FOC20).
Important: Ensure that the connectors are plugged in correctly (amplifier 1 into fiber optic 1, amplifier 2 into fiber optic 2.)



5


The amplifier is operated using the rechargeable battery, PowerPack.

Please charge the PowerPack before use! This may take up to 17 hours.

For details, please refer to the operating instructions of the amplifier which can be found on the Application Suite DVD.

6

Switch on your amplifier. The power LED must light up.



7

Connecting electrodes or electrode caps: Use the ribbon cable (BAC30 or BAC100) for connection.



8

After use: Disconnect the PowerPack from the amplifier and immediately recharge the PowerPack. Leave the PowerPack connected to the charger if not in use. The PowerPack can be trickle charged and overcharging is not possible.

Set Up Your BrainAmp System



Version 008 11/2015
© Brain Products GmbH

